



e-Brochure





Our company specializes in the manufacturing of ergonomically designed dental apparatus/furnishings /chairs etc. The company was established in year 2002 with a vision to become a global leader in our chosen field by endeavoring to become a recognized organization in the areas of durability, design, and quality.

For such vision to be accomplished, our valued and experienced management team has created an enjoyable and conducive working environment that stimulates, encourages, recognizes and respects each employee, instilling confidence and a positive attitude to all. As a result, it has produced an environment that develops and delivers the new features which are innovative and added improvements that are much better than our major competitors. As our highly competitive market is becoming ever-more design and quality focused, we can say confidently from our experience and from our employees' positive feedback that, we have the perfect winning formulas. The formula to attract best people of our choice and the formula to produce products that stands out as the best of the best.

2002. 8	- Establishment
2003. 7	- Introducing "Kaiser" in the market
2004. 1	- CE, ISO 9001, 13485 Certified.
2004. 3	- Introducing Dental unit & chair for Children "Le Petit Prince" in the market
2005. 3	- Recorded 1,000 set manufacturing of Kaiser



Kaiser

Kaiser, the most highly engineered dental unit & chair provides an advanced equipment range with maximal ergonomics and ease of use for the dentist.



Space art

Meet the most advanced chair in dental design and create a new space contouring the body.

With the highest standard of quality, design and value, your space will step up in time. This is the chair of the future....





Comfort in an ergonomic design

Visiting the dentist can be an uneasy experience for any patient. In order to ease the patient's fear, choose kaiser's advanced and safe technology. Experience the design with a comfortable and supportive headrest, ergonomical backrest and automatic safety switch.



Chair

HEAD REST

Double-articulating headrest cradle patient's head and neck. It also provides virtually limitless positioning for improved access when giving treatment even pedodontic treatment.

BACK REST

An ergonomically slim backrest is designed to provide the ultimate in access and comfort for the dentist's legs and knees.

ARM REST

A simple yet functional armrest provides easy access to the chair without hindrance. It gives optimal support and comfort for the patient during treatment.

SAFETY SWITCH

An automatic safety switch mounted on the bottom of the chair stops the chair from descending, preventing the risk of injury to the operator or any entrapped equipment.



Hygiene 365 days of the year

*With many patients visiting your clinic everyday,
hygiene is of the most importance.
The unit features a ceramic spittoon and suction filter,
both easily cleaned and maintained.
Kaiser's hygienic unit aids the dentist a relaxed,
healthy environment to work in.*



Unit

SPITTOON

A spittoon made of ceramic ensures superior protection against chipping and staining. Easily removable, the spittoon can be cleaned efficiently with hygiene in mind.

ASSISTANCE TABLE

With a smooth pivot point, the assistance table swings around the back of the unit, ensuring the assistant various access positions.

SUCTION FILTER

The suction filter is easily removed for cleaning and disinfection without needing to open the unit.





Simple design with a multifunctional display

A unit convenient to the dentist, this is the secret to efficient and easy work.

Treating patients can be long, tiring and harmful to the wrist.

Kaiser's instrument holder and multifunctional display control panel makes it the most user friendly unit available.



Cart Type

Table

INSTRUMENT HOLDER AND HANDPIECE TUBING

When selecting an instrument, other instruments are automatically turned off. The handpiece tubings are easily removed for cleaning and disinfection.

FILM VIEWER

The film viewer can be installed and removed effortlessly.

LCD CONTROL PANEL

A visual display unit on the keyboard shows the instrument in use and all selected functions.

The maximum and minimum light intensity can be selected by a push button. The LCD monitor then indicates the both light intensity modes. The built in microcomputer make three auto positions available for three user respectively. The preferred chair position can also be adjusted by an extra manual push button.



Table

INSTRUMENT HOLDER

*Instrument holder and table are fitted into one complete structure.
It is sturdily built yet has stylish modern appearance.*

FILM VIEWER

The film viewer can be installed and removed effortlessly without a power cable.

CONTROL PANEL

*The built in microcomputer make three auto positions available.
The preferred chair position can also be adjusted by an extra manual push button.
It is built to make the users utilize the ergonomically design touch board membrane.
Furthermore it is engineered to extend durability.
All the features have been put in place for the easy of use by the operator.*





Accessories



Dental light



Dr. Stool



LCD monitor




Distilled water tank



Assist Stool



Water warmer



Panorama film viewer



Dental film viewer



Foot controller



Piezo Scaler



Standard Accessories

3way syringe	1
Saliva Ejector	1
Operating Light	1

Options

- High speed handpiece
- Fiber optic handpiece
- Low speed handpiece
- Scaler
- Dental film viewer
- Panorama film viewer
- Doctor stool
- Water warmer
- Distilled water tank
- Curing light
- LCD monitor
- Suction(H.V.E)

Composition

- Motor type : Hydraulic Motor
- Chair Hight : 790° ±10 mm Max,
450° ±10 mm Min.
- Chair Length : 1,850 mm
- Light : 15,000 ~ 25,000 Lux

Kaiser Series

K700, K600, K500, K400

All the specification can be changed for improvement.

Color




Blue
(SD810)



Sky blue
(SD800)



Violet
(SD900)



Yellow
(SD200)



Pink
(SD500)



Green
(SD700)



111 AnyangMegavalley 799, Kwanyang-dong, Dongan-gu, Anyang, Kyungkido, Korea
TEL : (82 31)420-5710 FAX : (82 31)420-5987 MAIL: sky@skydental.co.kr

www.skydental.co.kr